

View Online at <https://aerobasegroup.com/nsn/6530-01-580-9432>

**Furnished Items:**

100 units; spore-g.Stearothermophilus; thermoplastic culture tube housing glass ampule and filter-paper carrier

**End Application:**

Test equipment for sterilization

**Protective Storage Feature:**

Controlled room temperature

**Package Model Number:**

Eztest h2o2 (for hydrogen peroxide)

**Usage Design:**

Sterilization biological indicator

**Approximate Shelf Life:**

24.0 months

**Product Name:**

Eztest h2o2 (for hydrogen peroxide; self-contained biological indicator

**Special Features:**

Bacterial spores on a filter-paper carrier packaged in a small, thrmoplastic culture tube; inside the tube is a sealed-glass ampoule of specially-formulated soybean casein digest culture medium containing a color indicator-turns yellow when spores grow

**Functional Description:**

Ten units per cycle is commonly used; expose load to sterilization cycle; with one non-sterilized unit mark positive control and the one exposed, crush glass culture ampoule to active; place in incubator 55-60 degrees c 24 hours; if all spores are killed the culture medium will remain red; if spores survived medium will turn yellow (positive)

**Special Test Features:**

For hydrogen peroxide-results in 24 hours without laboratory transfers

**Supplementary Features:**

Expiration is 24 months from date of manufacture; read results every 12 to 24 hours; all positive (yellow) samples should be recorded and discarded immediately. All yellow samples should be autoclaved prior to disposal

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fiig:**

A23800