

View Online at <https://aerobasegroup.com/nsn/4810-01-541-9847>

**Flow Control Device Normal Operation Position:**

Closed

**End Application:**

Gfe procurement lha/lhd classes

**Current Rating In Amps:**

5.000

**Valve Operation Method:**

Solenoid

**Voltage In Volts And Current Type:**

110.0 ac

**Product Name:**

Valve assembly, pinch; inavy: pinch valve

**Special Features:**

Next higher assembly: 6640-01-491-9330 auto blood glycer-/deglycerolization system (cell washing machine) model acp 215, used on machines sold beginning 2001

**Frequency In Hertz:**

50.0

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fig:**

A046b0