

View Online at <https://aerobasegroup.com/nsn/6505-01-414-9695>

**Usage Design:**

For the early detection of cancer consisting of the application of a special staining technique to cells exfoliated from diseased tissue

**Test Capability Quantity:**

1

**Component Breakdown Within Each Test Kit:**

Rpap smear tray, 1; disposable vaginal speculum, 1; 6 inch cotton tipped applicators, 2; 5 g sterile lubricant packet, 1; slide holder and mailer, 1; cervical scraper, 1; vinyl gloves, cuffed, 2; slides, frosted end, 2

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fiig:**

A510v0