

View Online at <https://aerobasegroup.com/nsn/6515-01-213-5989>

Special Features:

Co-set, closed injectate delivery system for cold injectate; used w/catheter intravenously to get cardiac output measurements by thermodilution; for use w/the american edwards 9520a cardiac output computer and a swanganz thermodilution catheter; after the system is primed and the cooling immersed in an ice bath, the injectate will reach operating temperature within 5 minutes; useable for up to 8 hours; 8 sets per case

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fig:

A23800