

View Online at <https://aerobasegroup.com/nsn/6515-01-230-9936>

Special Features:

Sterile; 25 gage; single use; formed; used for performing the anterior capsulotomy used w/5 ml syringe; may be attached to syringe filled w/air; balance salt solution to maintain the anterior chamber during capsulotomy; ; for opening a hole in condensed vitreous face; useful in dislocating the nucleus into the anterior chamber; bent to conform w/ the convexity of the crystalline lens; unity q type 5

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fig:

T133-c