

View Online at <https://aerobasegroup.com/nsn/6515-01-327-0717>

Usage Design:

Used on a short-term basis, percutaneously, in chest area to expand tissue in preparation for plastic surgery during breast reconstruction procedures

Material And Location:

Silicone

Size Designator:

600cc volume and 13cm width and 7cm projection and 20cm hemispheric dome

Features Provided:

Disposable

Special Features:

Sterile; with remote valve, round; can be used in right or left side of chest; color clear; dimensions are measurements of mandrels

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fig:

T133-a