

View Online at <https://aerobasegroup.com/nsn/6515-01-327-0717>

**Usage Design:**

Used on a short-term basis, percutaneously, in chest area to expand tissue in preparation for plastic surgery during breast reconstruction procedures

**Material And Location:**

Silicone

**Size Designator:**

600cc volume and 13cm width and 7cm projection and 20cm hemispheric dome

**Features Provided:**

Disposable

**Special Features:**

Sterile; with remote valve, round; can be used in right or left side of chest; color clear; dimensions are measurements of mandrels

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fig:**

T133-a