Hemodialysis Tubing Assembly - Page 1 of 1



View Online at https://aerobasegroup.com/nsn/6515-01-328-1858

Sterility:

Sterile

Disposition After Initial Use:

Disposable

Special Features:

Pre-pump arterial drip chamber; arterial and venous blood lines; protective caps; universal dialyzer connector (standard or locking);

priming cap; color coded locking patient connectors and slide clamps; protective injection sites; use with drake willock 480; 24 sets per case

Shelf Life:

N/a

Unit Of Measure:

Demilitarization:

No

Fiig:

A23800