

View Online at <https://aerobasegroup.com/nsn/6515-01-330-8171>

**Sterility:**

Sterile

**Disposition After Initial Use:**

Disposable

**Special Features:**

Arterial line consists of 2 injection sites, luer lock patient connector, luer lock dialyzer connector, 10 cm. Saline line, 40 cm. Heparin line, 30 cm. Nonitor line, pump segment; venous line consists of 2 injection sites, luer lock patient connector, luer lock patient connector, luer lock dialyzer connector, 30 cm. Monitor line; 124 ml. Total priming volume; for use with cobe centry ii, baxter sps, and drake-willock 480; 24 sets per case

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fiig:**

A23800