

View Online at <https://aerobasegroup.com/nsn/6515-01-330-8171>

Sterility:

Sterile

Disposition After Initial Use:

Disposable

Special Features:

Arterial line consists of 2 injection sites, luer lock patient connector, luer lock dialyzer connector, 10 cm. Saline line, 40 cm. Heparin line, 30 cm. Nonitor line, pump segment; venous line consists of 2 injection sites, luer lock patient connector, luer lock patient connector, luer lock dialyzer connector, 30 cm. Monitor line; 124 ml. Total priming volume; for use with cobe centry ii, baxter sps, and drake-willock 480; 24 sets per case

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

A23800