

View Online at <https://aerobasegroup.com/nsn/6515-01-369-5727>

Thickness:

5.0 millimeters

Overall Diameter:

15.0 millimeters

End Application:

Oes, evis and ultrasound scopes, olympus corp

Material And Location:

Rubber

Special Features:

Semi-disposable forceps valve no. 13; black; used as a seal/cap that fits over the biopsy valve port to prevent fluids from escaping; it may be punctured to extract fluids w/a syringe needle

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fig:

T133-h