

View Online at <https://aerobasegroup.com/nsn/6515-01-380-1650>

Disposition After Initial Use:

Disposable

Material And Location:

Plastic

Operation Method:

Plastic bulb

Nasal Tip:

Not provided

Special Features:

Assembled; adult mask with reinflation valve; low resistance duck bill valve; swivel patient connector reduces ett torque; star-lumen oxygen tubing with universal flow meter adapter; exhalation port; textured silicone-like bag performs in all climate conditions; for single patient use

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

T133-k