

View Online at <https://aerobasegroup.com/nsn/6515-01-380-1650>

**Disposition After Initial Use:**

Disposable

**Material And Location:**

Plastic

**Operation Method:**

Plastic bulb

**Nasal Tip:**

Not provided

**Special Features:**

Assembled; adult mask with reinflation valve; low resistance duck bill valve; swivel patient connector reduces ett torque; star-lumen oxygen tubing with universal flow meter adapter; exhalation port; textured silicone-like bag performs in all climate conditions; for single patient use

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fiig:**

T133-k