

View Online at <https://aerobasegroup.com/nsn/6515-01-447-4829>

**Attachment Method:**

Movable circumferential loop

**Usage Design:**

Elbow

**Size Designation:**

X-large adult

**Strap Quantity:**

1

**Usage Location:**

Either elbow

**Battened Section Fastening Method:**

Elastic strap w/snap fastener

**Special Features:**

W/tear shaped pad which can be placed for treatment of lateral or medial epicondylitis; pad is held down by a stretchable loop strap adjusting circumference to forearm; used for bursitis; carpal tunnel; contusions; osteoarthritis; sprain; tenosynovitis; bio-skin matl is lycra/fleece knit fabric

**Color:**

Black

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fiig:**

T133-j