NSN 6515-01-464-1670

Cardiovascular Catheterization Kit - Page 1 of 1



View Online at https://aerobasegroup.com/nsn/6515-01-464-1670

Nonsupply Items And Quantities:

Pmulti-lumen indwelling catheter: 12 french by 8 inches, radiopaque polyurethane with heparin coating, blue flextip, extension line clamps, injection site caps, 1; spring-wire guide, marked: 0.035 inch diameter by 26-13/16 inch, straight soft tip on one end - j tip on other, with arrow advancer, 1; needle: 18 gage by 2-1/2 inch, 5cc arrow roulerson spring-wire introduction syringe, 1; catheter: 18 gage by 2-1/2 inch radiopaque over 20 gage rw introducer needle, 1; pressure transduction probe: 20 gage by 3-1/8 inch tw, 1; needle: 25 gage by 1 inch, 3cc luer-lock syringe, 2; needle: 22 gage by 1-1/2 inch, 5cc luer-slip syringe, 1; 0.033 oz foil pack, 10 percent povidone-iodine ointment, 1; 3 swabstick foil pack, 10 percent povidone-iodine, 1; 5ml ampule hci, 1 percent lidocaine solution, 1; 24 by 36 inch drape with 4 inch fenestrations, 1; csr wrap, 1; 4 by 4 inch gauze pads, 5; 2 by 2 inch gauze pads, 2; 3-o silk suture, 1; no.11 scalpel, sharpsaway cup, dilator vessel, 1

inch, 3cc luer-lock syringe, 2; needle: 22 gage by 1-1/2 inch, 5cc luer-slip syringe, 1; 0.03
ointment, 1; 3 swabstick foil pack, 10 percent povidone-iodine, 1; 5ml ampule hci, 1 perce
with 4 inch fenestrations, 1; csr wrap, 1; 4 by 4 inch gauze pads, 5; 2 by 2 inch gauze pad
sharpsaway cup, dilator vessel, 1
Sterility:
Sterile
Disposition After Initial Use:
Disposable
Special Features:
Arrow-howes, large bore, multi-lumen central venous catheterization kit; 5 kits per case
Shelf Life:
N/a
Unit Of Measure:
Demilitarization:
No
Fiig:
A 22800