

View Online at <https://aerobasegroup.com/nsn/6515-01-491-2186>

**End Application:**

For use with the ent pocket light

**Product Name:**

Loop cerumen light child crs use w/ent pocket light 10s

**Special Features:**

Used to remove dirt, debris, or other small items from the ear canal of young children; with one hand the operator has a bright light beam focusing down the shaft into the ear canal leaving the other hand free

**Material:**

Stainless steel

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No