

View Online at <https://aerobasegroup.com/nsn/6515-01-492-6977>

**Usage Design:**

For plastic surgery procedures

**Jaw Grasping Method:**

Teeth

**Grasping Facility Quantity Per Jaw:**

1 first jaw and 2 second jaw

**Jaw Length:**

5.000 inches

**Special Features:**

Thumb tissue forceps

**Material And Location:**

Stainless steel

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fiig:**

A51710