

View Online at <https://aerobasegroup.com/nsn/6515-01-499-4076>

Handle Type:

Center action

Usage Design:

For peripheral iridectomies

Action Type:

Direct

Tip Type:

Angled and fine teeth

Unit Type:

Grasping

Grip Type:

Thumb-forefinger

Shank Type:

Curved

Jaw Type:

Angular

Jaw Grasping Method:

Fine teeth

Jaw Grasping Facility Location:

At tip

Shape:

Angular

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fig:

A51710