

View Online at <https://aerobasegroup.com/nsn/6515-01-503-3956>

Current Type:

Ac

End Application:

For use with the tps drill set; for use during mandibular/craniofacial surgery

Features Provided:

Reusable

Product Name:

Drill bone elect 5/32" ss jacobs chuck u/w tps drill set

Special Features:

5/32 inch drill; stainless steel; assists in accurate placement of drill bits and helps to prevent spooling of soft tissue

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

T133-e