

View Online at <https://aerobasegroup.com/nsn/6515-01-511-4706>

**Inside Diameter:**

0.4 millimeters

**Overall Length:**

6.0 millimeters

**Usage Design:**

Used as inner ear implant in stapedotomy procedures

**Special Features:**

Cause piston; stapes prosthesis; fluoroplastic; piston diameter 0.4 mm; loop inside diameter 0.6 mm; oval window niche; platinum armstrong ribbon

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fiig:**

T133-I