

View Online at <https://aerobasegroup.com/nsn/6515-01-521-4657>

Disposition After Initial Use:

Reusable

Product Name:

Connector saline lock/heparin lock hubs 100s

Special Features:

Saline lock (heparin lock hubs) will attach directly to any iv catheter hub; allowing needle safe access wiht any luer-lock syringe or luer-lock iv set; the connector is rated for 6 days of continuous use; swabable; latex free material; 100/pkg

Material:

Plastic

Shape:

Straight

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fig:

T117-a