

View Online at <https://aerobasegroup.com/nsn/6515-01-524-5800>

**Disposition After Initial Use:**

Reusable

**Product Name:**

Connector saline lock/heparin lock hub (1) per package

**Special Features:**

Saline lock (heparin lock hub) will attach directly to any iv catheter hub allowing needle safe access with any luer-lock syringe or luer-lock iv set; the connector is rated for 6 days of continuous use; swabable; latex free material

**Material:**

Plastic

**Shape:**

Straight

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fig:**

T117-a