

View Online at <https://aerobasegroup.com/nsn/6515-01-525-1499>

Usage Design:

Leg below knee

Size Designation:

Large

Material And Location:

Felt liner

Gender For Which Designed:

Male and female

Hardness Designation:

Firm

Strap Quantity:

5

Usage Location:

Either leg

Battened Section Fastening Method:

Web strap w/buckle

Extension Section Fastening Method:

Hook/loop closure

Special Features:

Motion control walker used following achilles tendon repairs, stable fractures of lower leg and ankle, acute ankle sprains or soft tissue injury; provides foot and ankle immobilization in 0, 10 or 20 degrees of plantarflexion as well as protected range of motion in 7.5 degree increments from 45 degrees plantarflexion to 30 degrees dorsiflexion

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fig:

T133-j