

View Online at <https://aerobasegroup.com/nsn/6515-01-573-0551>

Usage Design:

For use in the embolisation of brain arteriovenous malformations (avm's) with suitable angioarchitecture

Product Name:

Onyx avm system ld

Special Features:

System consists of 1.5 m. Vial of onyx 18 liquid embolic material, 1.5 ml vial of dms0 (dimethyl sulfoxide), three 1 ml delivery syringes

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

A23800