

View Online at <https://aerobasegroup.com/nsn/6515-01-573-5373>

**End Application:**

Used on edith flex heat and moisture exchanger for intubation anesthesia tracheotomy patients

**Special Features:**

Y filter/s plus edith flex; 30 mg oxygen moisture; removes 99.9 pct of bacteria/virus; 15f and 2 x 22m connectors; pkg of 20; in operating room, change after ea patient; in icu, change after 24 hrs

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fig:**

T117-a