Embolization Coil - Page 1 of 1



View Online at https://aerobasegroup.com/nsn/6515-01-587-4776

## End Application:

Used on artery/vein embolism surgery

## Material And Location:

Steel, corrosion resisting

Surface Treatment And Location:

Plastic, polytetrafluoroethylene

## **Special Features:**

Non stock special flipper detachable coil; do not use with polyvinyl or polyurethane catheters; recommend french size 5 or larger, 0.041 or

greater lumen, no sideports, nontapered multipurpose catheter

Shelf Life:

N/a

Unit Of Measure:

--

## Demilitarization:

No

Fiig:

T133-d