

View Online at <https://aerobasegroup.com/nsn/6515-01-657-8343>

Disposition After Initial Use:

Disposable

Product Name:

Airlife adapter valved tee 22 mm od/od, 30s

Special Features:

Used in conjunction with a nebulizer, the valved tee adapter offers a safe and convenient way to introduce medication into the circuit; its spring-loaded valve opens and closes the circuit during nebulizer insertion and removal; dim: 22mm i.D./o.D., latex free, single use, case of 30, case di (in) : l8.3 x h6.7

Shape:

T

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fig:

T117-a