

View Online at <https://aerobasegroup.com/nsn/6520-00-536-2770>

Physical Form:

Paste

Component Quantity:

2

Furnished Items And Quantity:

C/o one tube each of eugenol and zinc oxide

Setting Time:

Between 3 minutes and 8 minutes

Workable Consistency Temperature Rating:

68.0 degrees fahrenheit

Product Name:

Impression material dental zinc-oxide eugenol paste

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

T243-d