

View Online at <https://aerobasegroup.com/nsn/6520-01-062-7694>

**Product Name:**

Cement zinc oxide & eugenol dntl 15-20 min set time white

**Exterior Container Weight:**

6.00 pounds

**Unit Of Issue Cube:**

0.020 cubic feet

**General Description:**

100 grams; for use in lutting silver or gutta percha points; endodontically treated teeth when mixed with eugenol, 6505-00-153-8379; should have in vitro setting time of 15-20 minutes; in vivo working time should be at least one hour; grossman type 801

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No