

View Online at https://aerobasegroup.com/nsn/6520-01-062-7694

Product Name:

Cement zinc oxide & eugenol dntl 15-20 min set time white

Exterior Container Weight:

6.00 pounds

Unit Of Issue Cube:

0.020 cubic feet

General Description:

100 grams; for use in lutting silver or gutta percha points; endodontically treated teeth when mixed with eugenol, 6505-00-153-8379;

should have in vitro setting time of 15-20 minutes; in vivo working time should be at least one hour; grossman type 801

Shelf Life:

N/a

Unit Of Measure:

Demilitarization:

No