

View Online at <https://aerobasegroup.com/nsn/6520-01-333-0819>

**Usage Design:**

Model and denture purpose

**Setting Time:**

Between 3 minutes and 5 minutes

**Weight:**

25.0 pounds

**Product Name:**

Stone artificial dental 3-5min setting time 25lb white blue hard

**Special Features:**

White blue; hard; fast setting gypsum; for model anocclusal casts

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fiig:**

T243-d