## NSN 6520-01-377-2115

Dental Impression Material - Page 1 of 1



View Online at https://aerobasegroup.com/nsn/6520-01-377-2115

Physical Form:
Any acceptable
Component Quantity:
2
Setting Time:
4 minutes
Physical Form Color:
Purple
Workable Consistency Tempurature Rating:
75.0 degrees fahrenheit
Special Features:
1 tube base 80 ml; 1 tube catalyst 80 ml; mixing time 30 seconds; working time of 3 minutes; if refrigerated warm to room temperature
prior to use; delay model pouring for 3 hours after impression has taken
Material:
Plastic, vinyl polysiloxane
Material:
Plastic, vinyl polysiloxane
Shelf Life:
N/a
Unit Of Measure:
Demilitarization:
No
Fiig:
T243-d