

View Online at <https://aerobasegroup.com/nsn/6520-01-378-8901>

Physical Form:

Liquid and powder

Usage Design:

Dental restoration

Setting Time:

60 seconds

Curing Process Type:

Self

Special Features:

Amalgambond plus kit contains; 3 ml dentin activator (bottle "a"), 15 ml adhesive agent (bottle "aa"), 5 ml base (bottle "b"), 0.7 ml catalyst (syringe "c"); 2 gr high-performance additive (jar "hpa"); brushes, handles, tip , mixing well and foam pellets

Color:

Any acceptable

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

T243-d