

View Online at <https://aerobasegroup.com/nsn/6520-01-388-0445>

Physical Form:

Powder

Component Quantity:

1

Setting Time:

3 minutes

Weight:

4.500 pounds

Special Features:

A mixture of high grade plaster and natural gums; when mixed with water forms a paste with putty-like consistency; may be used with all edentulous techniques for full dentures, corrective impressions, and for rebasing dentures; may be used over compound, base plate trays, fastray, acrylic, and all new or old dentures; may be used for the mcgrane non-pressure technique or where pressure of non-pressure is indicated

Material:

Any acceptable

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

T243-d