

View Online at <https://aerobasegroup.com/nsn/6520-01-388-8597>

Usage Design:

Bridge restoration and crown restoration

Setting Time:

12 minutes

Product Name:

Investment dental bridge/crown restoration 12min setting time24s

Special Features:

Complete porcelain to metal, crown and bridge investment; working time 6 to 8 minutes; unit pkg c/o measuring vial, ring liner, instructions, 8 oz liquid and 24-90 gm pouches of powder

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

T243-d