Dental Impression Material - Page 1 of 1



View Online at https://aerobasegroup.com/nsn/6520-01-390-0747

Physical Form:

Powder

Component Quantity:

1

Setting Time:

Between 3 minutes and 6 minutes

Workable Consistency Tempurature Rating:

Between 68.0 degrees fahrenheit and 77.0 degrees fahrenheit

Hydrocolloid Base Type:

Alginate

Special Features:

Regular set; not dustless; 12 foil pouches/pg (21g each pouch)

Material:

Any acceptable

Shelf Life:

N/a

Unit Of Measure:

--Demilitarization:

Deminiarizat

No

Fiig:

T243-d