

View Online at <https://aerobasegroup.com/nsn/6520-01-435-1564>

**Operating Temperature Range:**

+0002.0/+1200.0 degrees fahrenheit

**Voltage Rating:**

Between 100.0 volts and 115.0 volts

**Product Name:**

Furnace dental laboratory elec used w/vacuum pump p/n 701794

**Special Features:**

Cphoenix quick cool porcelain furnace; firing temperature range 199 degrees to 2192 degrees f; an oilless vacuum pump draws hot air away from muffle quickly; microprocessor controlled vacuum start and stop

**Supplementary Features:**

200/240 v, 60hz available; microprocessor update capability;

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fig:**

T243-b