

View Online at <https://aerobasegroup.com/nsn/6520-01-435-1564>

Operating Temperature Range:

+0002.0/+1200.0 degrees fahrenheit

Voltage Rating:

Between 100.0 volts and 115.0 volts

Product Name:

Furnace dental laboratory elec used w/vacuum pump p/n 701794

Special Features:

Cphoenix quick cool porcelain furnace; firing temperature range 199 degrees to 2192 degrees f; an oilless vacuum pump draws hot air away from muffle quickly; microprocessor controlled vacuum start and stop

Supplementary Features:

200/240 v, 60hz available; microrocessor update capability;

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

T243-b