

View Online at <https://aerobasegroup.com/nsn/6520-01-468-4747>

Usage Design:

Orthodontist procedures

Special Features:

Nickel titanium wire, reverse curve of spee arch (rcs); recommended for bite opening/closing and to correct rcs irregularities; counterbalance for elastomeric chain; acts 3 dimensionally to forstall unwanted tooth movement; 0.016 in by 0.022 in, upper

Surface Treatment:

Nickel

Surface Treatment:

Nickel

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

T243-d