

View Online at <https://aerobasegroup.com/nsn/6520-01-468-4747>

**Usage Design:**

Orthodontist procedures

**Special Features:**

Nickel titanium wire, reverse curve of spee arch (rcs); recommended for bite opening/closing and to correct rcs irregularities; counterbalance for elastomeric chain; acts 3 dimensionally to forstall unwanted tooth movement; 0.016 in by 0.022 in, upper

**Surface Treatment:**

Nickel

**Surface Treatment:**

Nickel

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fiig:**

T243-d