

View Online at <https://aerobasegroup.com/nsn/6520-01-478-5899>

End Application:

Used on model km110

Err Default:

Err default

Special Features:

Contents: acetone and isobutyl acetate; flammable, avoid contact w/eyes and skin; may be applied to columbia dentoform teeth w/sm syring, toothpick or like applicator; matl dries in 3 to 4 hrs at room temp, fully cured in 24 hrs; 1 tube does about 24 models

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

T243-a