

View Online at <https://aerobasegroup.com/nsn/6520-01-533-3779>

**Product Name:**

Tip ultrasonic f/root canal cleaning sz buc number 2a

**Special Features:**

Disk-like radiused tip, can be used to smooth and safely plane attached pulp stones from the pulp chamber floor without scoring it; small version: 1.0mm diameter used for corners of molar access prep and bicuspid access prep; use at power level 1-5; size: buc-2a

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fig:**

T243-a