

View Online at <https://aerobasegroup.com/nsn/6520-01-543-8198>

Usage Design:

Orthodontist procedures

Curing Process Type:

Light

Special Features:

Mini mold starter kit, 3 tips and 2 handles, instructions; fabricate buttons, bite ramps; habit control spikes, temp brackets, molar tubes, pads and lingual retainer bonds using light cure adhesives; ea mold packs w/light cure adhesive; press against tooth, cure and peel away mold leaving formed auxiliary on any tooth in minutes; compatible w/quick cure adhesives and reliance assure light bond

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

T243-d