Denture Base-orthodontic Material - Page 1 of 1



View Online at https://aerobasegroup.com/nsn/6520-01-543-8198

#### **Usage Design:**

Orthodontist procedures

## **Curing Process Type:**

Light

### **Special Features:**

Mini mold starter kit, 3 tips and 2 handles, instructions; fabricate buttons, bite ramps; habit control spikes, temp brackets, molar tubes, pads and lingual retainer bonds using light cure adhesives; ea mold packs w/light cure adhesive; press against tooth, cure and peel away mold leaving formed auxiliary on any tooth in minutes; compatible w/quick cure adhesives and reliance assure light bond **Shelf Life:** 

N/a

### Unit Of Measure:

--

# Demilitarization:

No

Fiig:

T243-d