

View Online at <https://aerobasegroup.com/nsn/6525-01-514-9962>

Mounting Type:

Stationary

Universal Medical Device Nomenclature System Code:

18429

Operation Method:

Automatic

Tube Quantity:

1

Tube Type:

Double focus, rotating anode radiographic

Table Type:

Nontilting

Rectification Type:

Full-wave

Tube Cooling Method:

Air

Tube Mounting Type:

Tube stand

Radiographic Timer Type:

Exposure

Transformer Rectification Method:

Self-rectifying

Product Name:

Bucky diagnostic x-ray unit

Phase:

Three

Special Features:

190 to 240vac; stationary free standing unit; solid state circuitry; rectified high voltage transformer and cables; free standing control cabinet with visual indicator; for field military medical facilities; includes: bucky unit, optimus 50 generator, patient table, wall mount operator panel, tube assembly panel, arm rest, bucky ve/vt, lateral cassette holder, bucky diagnostic ve/vt, high speed rotor control, automatic collimation

Frequency In Hertz:

50.0 or 60.0

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

T140-c