

View Online at <https://aerobasegroup.com/nsn/6530-00-079-7039>

**Physical Form:**

Strip

**Product Name:**

Indicator sterilization indicator color chg by steam 250s

**Special Features:**

Indicator color change by steam heat; center line changes color from white to brown when proper sterilization conditions of steam; temperature and time have occurred within pack, at 250 degf aafter 12 to 15 min. Exposure, and 270 degrees f after 2 to 3 min exposure; steam; item w/potency period or expiration date; unit qty, 250 strips

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fiig:**

T396-a