

View Online at <https://aerobasegroup.com/nsn/6530-01-213-9503>

**Special Features:**

Steam monitor strip can be used in both gravity/high vacuum sterilizers, at temperatures between 250-275 degree f; autoclaved; bar turns black to indicate exposure to pressurized steam; helps detect gasket leaks, trap leaks, vacuum malfunction, time or temperature control malfunctions, improper loading/improper wrapping; unit q type 250

**General Description:**

Steam monitor strip can be used in both gravity/high vacuum sterilizers, at temperatures between 250-275 degree f; autoclaved; bar turns black to indicate exposure to pressurized steam; helps detect gasket leaks, trap leaks, vacuum malfunction, time or temperature control malfunctions, improper loading/improper wrapping; unit q type 250

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fiig:**

T396-a