

View Online at <https://aerobasegroup.com/nsn/6530-01-224-6058>

Product Name:

Indicator steriliz disp bars turn blk indic satisfactory cycle 50

Special Features:

Disposable; monitoring bars on the dart indicator will turn uniformly black meaning cycle is satisfactory; uneven color change is an indication of failure; test results are clear cut; elimination the chance of misinterpretation + four minutes exposure time at 270 degrees f.

During a pre-vacuum sterilization cycle

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

T396-a