

View Online at <https://aerobasegroup.com/nsn/6530-01-224-6058>

**Product Name:**

Indicator steriliz disp bars turn blk indic satisfactory cycle 50

**Special Features:**

Disposable; monitoring bars on the dart indicator will turn uniformly black meaning cycle is satisfactory; uneven color change is an indication of failure; test results are clear cut; elimination the chance of misinterpretation + four minutes exposure time at 270 degrees f.

During a pre-vacuum sterilization cycle

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fiig:**

T396-a