

View Online at <https://aerobasegroup.com/nsn/6530-01-348-6068>

**Special Features:**

Dart designed to test for the removal of air from the chamber during the pre-vacuum stage; disposable; monitoring bars on the dart indicator will turn uniformly black meaning cycle is satisfactory; uneven color change is an indication of failure; test results are clear cut; eliminating the chance of misinterpretation four minutes exposure time at 270 degrees f during a pre-vacuum sterilization cycle

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fig:**

T396-a